

Product datasheet for **RC203434A1V**

Human PNO1 (NM_020143) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human PNO1 (NM_020143) AAV Particle
Tag:	Myc-DDK
Symbol:	PNO1
Synonyms:	KHRBP1; RRP20
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203434 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAATCCGAAATGGAACGCAGAGCGCCAGGGCAGAGGAGGGCTTTACCCAGGTACCCGCAAGGGTG
GCCGACGGGCGAAGAAACGACAGGCTGAACAGCTGTCCGACGAGGAGGGCGGGGATGCGGGCCGCAT
GGACACAGAGGAGGCCAGGCCGCGCAAGAGGCCCGTCTCCACCCCTCTGTGGGACGGGCTCCTGAGT
GGGAAAGAAGAAACAAGGAAATTCAGTCCAGCTAACAGATACACACCATTGAAAGAAAACCTGGATGA
AGATATTTACTCCTATTGTGGAACATTTGGGACTTCAGATACGCTTTAACTTGAAATCAAGGAATGTAGA
AATCAGGACTTGTAAAGAAACCAAGGATGTTAGTGCTCTGACAAAAGCAGCTGATTTTGTGAAAGCTTTT
ATTCTCGGCTTTTCAAGTGGAGGATGCACTTGCCCTCATCAGGTTGGATGACCTCTTCTAGAGTCTTTTG
AAATTACAGATGTTAAACCCCTAAAGGGAGACCATCTATCCAGGGCAATAGGAAGAATCGCTGGCAAAGG
AGGAAAAACCAATTCACCATAGAGAATGTGACACGGACAAGGATAGTTTGGCTGATGTGAAAGTTCAC
ATCCTTGGCTCCTTCCAAATATCAAGATGGCAAGAACTGCCATTTGCAACCTAATCTTGGGAAATCCTC
CTTCAAGGTTTATGGCAATATTCGAGCTGTGGCTAGCAGATCAGCAGATCGATTC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC203434 protein sequence Red=Cloning site Green=Tags(s) MESEMETQSARAEEGFTQVTRKGGRRAKKRQAEQLSAAGEGGDAGRMDTEEARPAKRPVFPPLCGDGLLS GKEETRKIPVPANRYTPLKENWMKIFTPIVEHLGLQIRFNLKSRNVEIRTCKETKDVSAITKAADFVKAF ILGFQVEDALALIRLDDLFLSFETIDVKPLKGDHLSRAIGRIAGKGGTKFTIENVTRTRIVLADVVKVH ILGSFQNIKMARTAICNLILGNPPSKVYGNIRAVASRSADRF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_020143
ORF Size:	756 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_020143.2</u>
RefSeq Size:	1994 bp
RefSeq ORF:	759 bp
Locus ID:	56902
UniProt ID:	<u>Q9NRX1</u>
Cytogenetics:	2p14
MW:	27.9 kDa