

Product datasheet for **RC203447L4V**

NGLY1 (NM_018297) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NGLY1 (NM_018297) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NGLY1
Synonyms:	CDDG; CDG1V; PNG-1; PNG1; PNGase
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018297
ORF Size:	1965 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203447).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018297.2
RefSeq Size:	2534 bp
RefSeq ORF:	1965 bp
Locus ID:	55768
UniProt ID:	Q96IV0
Cytogenetics:	3p24.2
Domains:	TGc, PUG, PAW
MW:	74.4 kDa



[View online »](#)

Gene Summary:

This gene encodes an enzyme that catalyzes hydrolysis of an N(4)-(acetyl-beta-D-glucosaminy) asparagine residue to N-acetyl-beta-D-glucosaminylamine and a peptide containing an aspartate residue. The encoded enzyme may play a role in the proteasome-mediated degradation of misfolded glycoproteins. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Feb 2009]