

Product datasheet for **RC203455A1V**

Human PICK1 (NM_001039584) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human PICK1 (NM_001039584) AAV Particle
Tag:	Myc-DDK
Symbol:	PICK1
Synonyms:	PICK; PRKCABP
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203455 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTTTGCAGACTTGGATTATGACATCGAAGAGGATAAACTCGGAATCCCGACTGTGCCTGGGAAGGTGA
CCCTGCAGAAGGATGCTCAGAACCTGATCGGGATCAGCATTGGAGGAGGGGCCAGTACTGTCCCTGCCT
CTATATCGTCCAGGTATTTGACAACACCCAGCAGCCTTGGACGGCACAGTGGCAGCTGGCGATGAGATC
ACCGGTGTCAATGGCAGGTCAATCAAAGGGAAAACTAAGGTGGAGGTGGCGAAGATGATTCAAGAGGTGA
AGGGGGAGGTGACCATCACTACAACAAGCTGCAGGCGGACCCCAAGCAGGGCATGTCCTGGACATTGT
GTTGAAGAAAGTCAAGCACCGGCTGGTGGAGAATGAGTTCAAGGACCGCAGATGCTCTGGGCCTGAGC
CGGGCCATCCTGTGCAATGATGGGCTTGTCAAGAGGCTAGAGGAGCTGGAGCGGACCGCTGAGCTATACA
AAGGGATGACGGAACACACCAAGAACCTCCTACGGGCCTTTTATGAGCTGTGCGAGACTACCGGGCCTT
TGGGGACGTGTTCTCCGTGATCGGGGTGCGGGAGCCCCAGCCAGCTGCGAGCGAGGCTTTTGTGAAGTTC
GCCGATGCCACCGCAGCATCGAGAAGTTCGGCATTGGCTTCTGAAAACCATCAAGCCGATGCTGACGG
ATCTGAACACGTACCTCAACAAAGCCATCCCGGACACTCGCCTCACCATCAAGAAGTACCTGGACGTGAA
GTTTGAGTACCTGTCGTAAGGTGAAGGAGATGGATGACGAGGAATACAGCTGCATTGCCCTA
GGCGAGCCCTTTACCGGTGAGCACCGCAACTATGAGTACCGCCTGATCCTGCGCTGCCGCCAGGAGG
CGCGCGCCCGCTTCTCCAGATGCGCAAGGATGTGCTGGAGAAGATGGAGCTGCTGGACCAGAAGCACGT
CCAGGACATCGTGTTCCAGCTGCAGCGCCTCGTGTCCACCATGTCCAAGTACTACAACGACTGCTACGCA
GTGCTGCGGGATGCCGACGTCTTCCCCATCGAGGTAGACCTGGCGCACACCACATTGGCCTATGGCCTCA
ACCAGGAGGAGTTCACAGATGGGGAGGAGGAGGAGGAGGAAGACACGGCAGCTGGGGAGCCGTCCAG
GGATACACGAGGGGCTGCTGGGCCCTTGACAAGGTTGAAGCTGGTGTGACTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC203455 protein sequence Red=Cloning site Green=Tags(s) MFADLDYDIEEDKLGIPVPGKVTLQKDAQNLIGISIGGGAQYCPCLYIVQVFDNTPAALDGTVAAGDEI TGVNGRSIKGKTKVEVAKMIQEVKGEVTIHYNKLQADPKQGMSLDIVLKKVKHRLVENMSSGTADALGLS RAILCNDGLVKRLEELERTAELYKGMTEHTKNLLRAFYELSQTHRAFGDVFSVIGVREPQPAASEAFVKF ADAHRSIEKFGIRLLKTIKMLTDLNTYLNKAIPDTRLTIKKYLDVKFEYLSYCLKVKEMDDEEYSCIAL GEPLYRVSTGNYEYRLILRCRQEARARFSQMRKDVLEKMELLDQKHVQDIVFQLQRLVSTMSKYND CYA VLRDADVPFPIEVDLAHTTLAYGLNQEEFTDGEEEEEEDTAAGEPSRDTRGAAGPLDKGGSWCDS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001039584
ORF Size:	1245 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001039584.1</u> , <u>NP_001034673.1</u>
RefSeq Size:	2073 bp
RefSeq ORF:	1248 bp
Locus ID:	9463
UniProt ID:	<u>Q9NRD5</u>
Cytogenetics:	22q13.1
MW:	46.6 kDa