

Product datasheet for **RC203462L3V**

POLG2 (NM_007215) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	POLG2 (NM_007215) Human Tagged ORF Clone Lentiviral Particle
Symbol:	POLG2
Synonyms:	HP55; MTDPS16; MTPOLB; PEOA4; POLB; POLG-BETA; POLGB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007215
ORF Size:	1455 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203462).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007215.2
RefSeq Size:	1607 bp
RefSeq ORF:	1458 bp
Locus ID:	11232
UniProt ID:	Q9UHN1
Cytogenetics:	17q23.3
Domains:	HGTP_anticodon
Protein Families:	Stem cell - Pluripotency



[View online »](#)

Protein Pathways: Metabolic pathways

MW: 54.9 kDa

Gene Summary: This gene encodes the processivity subunit of the mitochondrial DNA polymerase gamma. The encoded protein forms a heterotrimer containing one catalytic subunit and two processivity subunits. This protein enhances DNA binding and promotes processive DNA synthesis. Mutations in this gene result in autosomal dominant progressive external ophthalmoplegia with mitochondrial DNA deletions.[provided by RefSeq, Sep 2009]