

### Product datasheet for RC203463L1

# CUL4A (NM\_003589) Human Tagged Lenti ORF Clone

#### **Product data:**

**Product Type:** Expression Plasmids

Product Name: CUL4A (NM 003589) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: CUL4A

Mammalian Cell

Selection:

None

Vector: pLenti-C-Myc-DDK (PS100064)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clor


Sequence:

The ORF insert of this clone is exactly the same as(RC203463).

Restriction Sites: Sgfl-Mlul

**Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF.

**ACCN:** NM\_003589 **ORF Size:** 2280 bp



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



#### CUL4A (NM\_003589) Human Tagged Lenti ORF Clone - RC203463L1

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

**RefSeq:** <u>NM 003589.2</u>

 RefSeq Size:
 3735 bp

 RefSeq ORF:
 1980 bp

 Locus ID:
 8451

UniProt ID: Q13619
Cytogenetics: 13q34

Domains: CULLIN

Protein Pathways: Nucleotide excision repair, Ubiquitin mediated proteolysis

**MW:** 87.7 kDa


**Gene Summary:** CUL4A is the ubiquitin ligase component of a multimeric complex involved in the degradation

of DNA damage-response proteins (Liu et al., 2009 [PubMed 19481525]).[supplied by OMIM,


Oct 20091



## **Product images:**



Circular map for RC203463L1



Double digestion of RC203463L1 using Sgfl and Mlul  $\,$