

Product datasheet for **RC203469L4V**

ARL3 (NM_004311) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ARL3 (NM_004311) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ARL3
Synonyms:	ARFL3; JBTS35; RP83
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004311
ORF Size:	546 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203469).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004311.2
RefSeq Size:	3891 bp
RefSeq ORF:	549 bp
Locus ID:	403
UniProt ID:	P36405
Cytogenetics:	10q24.32
Domains:	SAR, ARF, arf
Protein Families:	Druggable Genome



[View online »](#)

MW: 20.5 kDa

Gene Summary: ADP-ribosylation factor-like 3 is a member of the ADP-ribosylation factor family of GTP-binding proteins. ARL3 binds guanine nucleotides but lacks ADP-ribosylation factor activity. [provided by RefSeq, Jul 2008]