

Product datasheet for **RC203477A1V**

Human TFIIE beta (GTF2E2) (NM_002095) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TFIIE beta (GTF2E2) (NM_002095) AAV Particle
Tag:	Myc-DDK
Symbol:	TFIIE beta
Synonyms:	FE; TF2E2; TFIIE-B; TTD6
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203477 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGATCCAAGCCTGTTGAGAGAAAGGGAGCTGTTCAAAAACGAGCTCTTTCTACTCCTGTAGTAGAAA
AACGTTTCAGCATCTTCTGAGTCATCATCATCATCGTCAAAGAAGAAGAAAACAAAGGTAGAACATGGAGG
ATCGTCAGGCTCTAAACAAAATTCTGATCATAGCAATGGATCATTTAACTTGAAGCTTTGTCAGGAAGC
TCTGGATATAAGTTTGGTGTCTTGCTAAGATTGTGAATTACATGAAGACACGGCATCAGCGAGGAGATA
CGCATCCTCTAACCTTAGATGAAATTTTGGATGAAACACAACATTTAGATATTGGACTCAAGCAGAAACA
ATGGCTAATGACTGAGGCTTTAGTCAACAATCCCAAAATTGAAGTAATAGATGGGAAGTATGCTTTCAAG
CCCAAGTACAACGTGAGAGATAAGAAGGCCCTACTTAGGCTCTTAGATCAGCATGACCAGCGAGGATTAG
GAGGAATTCCTTTAGAAGACATAGAAGAAGCACTGCCCAATCCAGAAAGCTGTCAAGGCTTTGGGGGA
CCAGATACTATTTGTAATCGTCCCGATAAGAAGAAAATACTTTCTTCAATGATAAGAGCTGTCAGTTT
TCTGTGGATGAAGAATTCAGAACTGTGGAGGAGTGCACTGTAGATTCATGGACGAGGAGAGAAATG
AAGAATATCTGAAGCGACAGGATTTTCTTCCATGCAGGAATCTGGACCAAGAAAGTGGCCCTATTCA
GAGAAGGAAAAAGCCTGCTTACAGAAAAAGCGACGCTTAAGACTCATAACGAACACTTGGCTGGAGTG
CTGAAGGATTACTCTGACATTACTTCCAGCAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC203477 protein sequence Red=Cloning site Green=Tags(s) MDPSLLRERELFKKRALSTPVVEKRSASSESSSSSSSKKKKTKVEHGGSSGSKQNSDHSNGSFNLKALSGS SGYKFGVLAKIVNYMKTRHQRGDTHPLTLDEILDETQHLDIGLKQKQWLMTEALVNNPKIEVIDGKYAFK PKYNVRDKKALLRLDQHDQRLGGILLEDIEEALPNSQKAVKALGDQILFVNRPDKKKILFFNDKSCQF SVDEEFQKLWRSVTVDSDMEEKIEEYLKRQGISSMQESGPKKVAPIQRRKKPASQKKRRFKTHNEHLAGV LKDYSIDITSSK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_002095
ORF Size:	873 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_002095.4</u>
RefSeq Size:	1560 bp
RefSeq ORF:	876 bp
Locus ID:	2961
UniProt ID:	<u>P29084</u>
Cytogenetics:	8p12
MW:	33 kDa