

Product datasheet for **RC203503L4V**

KRTHA3B (KRT33B) (NM_002279) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KRTHA3B (KRT33B) (NM_002279) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KRTHA3B
Synonyms:	Ha-3II; HA3II; hHa3-II; K33B; KRTHA3A; KRTHA3B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002279
ORF Size:	1212 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203503).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002279.3
RefSeq Size:	1626 bp
RefSeq ORF:	1215 bp
Locus ID:	3884
UniProt ID:	Q14525
Cytogenetics:	17q21.2
MW:	46.2 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the keratin gene family. This gene is one of multiple type I hair keratin genes that are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription. As a type I hair keratin, the encoded protein is an acidic protein which heterodimerizes with type II keratins to form hair and nails. There are two isoforms of this protein, encoded by two separate genes, keratin 33A and keratin 33B. [provided by RefSeq, May 2012]