

Product datasheet for **RC203505L3V**

Clathrin light chain (CLTA) (NM_001833) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Clathrin light chain (CLTA) (NM_001833) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Clathrin light chain
Synonyms:	LCA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001833
ORF Size:	654 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203505).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001833.2
RefSeq Size:	1156 bp
RefSeq ORF:	657 bp
Locus ID:	1211
UniProt ID:	P09496
Cytogenetics:	9p13.3
Domains:	Clathrin_Ig_ch
Protein Pathways:	Endocytosis, Huntington's disease, Lysosome



[View online »](#)

MW: 23.7 kDa

Gene Summary: Clathrin is a large, soluble protein composed of heavy and light chains. It functions as the main structural component of the lattice-type cytoplasmic face of coated pits and vesicles which entrap specific macromolecules during receptor-mediated endocytosis. This gene encodes one of two clathrin light chain proteins which are believed to function as regulatory elements. Alternative splicing results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 8 and 12. [provided by RefSeq, May 2010]