

Product datasheet for **RC203524L4V**

TCEAL3 (NM_001006933) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | TCEAL3 (NM_001006933) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | TCEAL3 |
| Synonyms: | WEX8 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001006933 |
| ORF Size: | 600 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC203524). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001006933.1 , NP_001006934.1 |
| RefSeq Size: | 1154 bp |
| RefSeq ORF: | 603 bp |
| Locus ID: | 85012 |
| UniProt ID: | Q969E4 |
| Cytogenetics: | Xq22.2 |
| MW: | 22.5 kDa |


[View online »](#)

Gene Summary:

This gene encodes a member of the transcription elongation factor A (SII)-like (TCEAL) gene family. Members of this family contain TFA domains and may function as nuclear phosphoproteins that modulate transcription in a promoter context-dependent manner. Multiple family members are located on the X chromosome. Alternative splicing results in multiple transcript variants encoding a single isoform. [provided by RefSeq, Jul 2008]