

Product datasheet for **RC203525L4V**

DACT3 (NM_145056) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DACT3 (NM_145056) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DACT3
Synonyms:	DAPPER3; RRR1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145056
ORF Size:	1212 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203525).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145056.1 , NP_659493.1
RefSeq Size:	2834 bp
RefSeq ORF:	1890 bp
Locus ID:	147906
UniProt ID:	Q96B18
Cytogenetics:	19q13.32
MW:	41.5 kDa



[View online »](#)

Gene Summary:

May be involved in regulation of intracellular signaling pathways during development. Specifically thought to play a role in canonical and/or non-canonical Wnt signaling pathways through interaction with DSH (Dishevelled) family proteins.[UniProtKB/Swiss-Prot Function]