

## Product datasheet for **RC203543L1V**

### NUDC (NM\_006600) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	NUDC (NM_006600) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NUDC
Synonyms:	HNUDC; MNUDC; NPD011
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_006600
ORF Size:	993 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203543).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006600.2</a>
RefSeq Size:	1813 bp
RefSeq ORF:	996 bp
Locus ID:	10726
UniProt ID:	<a href="#">Q9Y266</a>
Cytogenetics:	1p36.11
Domains:	Nuclear_move
MW:	38.2 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a nuclear distribution protein that plays an essential role in mitosis and cytokinesis. The encoded protein is involved in spindle formation during mitosis and in microtubule organization during cytokinesis. Pseudogenes of this gene are found on chromosome 2. [provided by RefSeq, Feb 2012]