

Product datasheet for **RC203547L4V**

ZNF213 (NM_004220) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF213 (NM_004220) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF213
Synonyms:	CR53; ZKSCAN21; ZSCAN53
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004220
ORF Size:	1377 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203547).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004220.1
RefSeq Size:	3301 bp
RefSeq ORF:	1380 bp
Locus ID:	7760
UniProt ID:	O14771
Cytogenetics:	16p13.3
Protein Families:	Transcription Factors
MW:	51.3 kDa



[View online »](#)

Gene Summary:

C2H2 zinc finger proteins, such as ZNF213, have bipartite structures in which one domain binds DNA or RNA and the other modulates target gene expression.[supplied by OMIM, Apr 2004]