

Product datasheet for **RC203550L3V**

HPDL (NM_032756) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HPDL (NM_032756) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HPDL
Synonyms:	4-HPPD-L; GLOXD1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032756
ORF Size:	1113 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203550).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032756.2
RefSeq Size:	1803 bp
RefSeq ORF:	1116 bp
Locus ID:	84842
UniProt ID:	Q96IR7
Cytogenetics:	1p34.1
MW:	39.4 kDa


[View online »](#)

Gene Summary:

The protein encoded by this intronless gene localizes to mitochondria, where it may function as 4-hydroxyphenylpyruvate dioxygenase. Clinical studies have identified several bi-allelic variants in this gene that lower the level of the encoded protein and lead to a clinically variable form of pediatric-onset spastic movement disorder. [provided by RefSeq, Aug 2020]