

Product datasheet for **RC203563A1V**

Human TBCC (NM_003192) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TBCC (NM_003192) AAV Particle
Tag:	Myc-DDK
Symbol:	TBCC
Synonyms:	CFC
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203563 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGTCCGTCACTTGGCTCCGCTGCTGCTCAGGACCGGAGACATGGAGTCCCAGCGGGACCTGAGCC
TGGTGCCTGAGCGGCTTCAGAGACGCGAACAAGAACGGCAGCTGGAAGTTGAAAGGCGGAAACAAAAGCG
GCAGAACAGGAGGTAGAGAAGGAGAACAGCCACTTTTCGTCGCCACCTTTGCTCGGGAGCGAGCGGCC
GTGGAAGAGCTTCTGGAGCGCGGAGTCGGTCGAGCGGCTGGAGGAGCGGCCCTCTCGGCTCCAGGGGC
TGCAGAACTAATCAACGACTCAGTTTTTTTCTAGCCGCTTACGACCTGCGGCAGGGACAAGAGGCGCT
GGCGCGGCTGCAGGCGGCCCTTGCCGAGCGCGCGCGGGGCTGCAGCCCAAGAAGCGTTTCGCTTTCAAG
ACCCGGGGAAGGATGCTGCTTCTACCAAGTAGACGCGGCTCCTGGCATCCCCCGGCAGTTGAAA
GCATACAGGACTCCCCGCTGCCAAGAAGGCGGAAGGAGACCTCGGCCCGAGCTGGGTCTGCGGTTTCTC
CAACCTGGAGTCCCAAGTCTTGAGAAGAGAGCCAGCGAGTTGCACCAGCGGACGTTCTTTGACCGAA
CTGAGCAACTGCACGGTCAGACTTTATGGAATCCCAACACCCTGCGGCTAACCAAGGCCACAGCTGCA
AGCTGCTCTGCGGTCCGGTGTCTACCTCTGTTTCTTGAGGAGTGCAGTGACTGCGTGTGGCAGTGGC
CTGCCAACAGCTCCGCATACACAGTACGAAAGACACCCGATCTTCTGCAGGTGACCAGCAGGGCCATC
GTGGAGGACTGCAGTGGGATCCAGTTTCCCCCTTACACCTGGAGTACCCGAGATCGACAAGGACTTCG
AGAGCTCTGTTTAGATAGGAGCAAAAATAACTGGAACGATGTTGACGATTTTAACTGGCTGGCCCGGA
TATGGCTCCCCAACTGGAGTATTCTTCTGAAGAGGAGCGAAATATCCAGTGGGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC203563 protein sequence Red=Cloning site Green=Tags(s) MESVSCSAAAVRTGDMESQRDL SLVPERLQRREQERQLEVERRKQKRQNEVEKENS HFFVATFARERAA VEELLERAESVERLEEAASRLQGLQKL INDSVFFLAAYDLRQGQEQALARLQAALAERRRGLQPKKRFAFK TRGKDAASSTKVDAAPGIPPAVESIQDSPLPKKAEGDLGPSWVCGF SNLESQVLEKRASELHQRDVLLTE LSNCTVRLYGNPNTLR LTKAHSCKLLCGPVSTSVFLEDCSDCVLAVACQQLRIHSTKDTRIFLQVTSRAI VEDCSGIQFAPYTWSYPEIDKDFESSGLDRSKNNWNVDDFNWLARDMASPNWSILPEEERNIQWD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_003192
ORF Size:	1038 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_003192.2</u>
RefSeq Size:	1675 bp
RefSeq ORF:	1041 bp
Locus ID:	6903
UniProt ID:	<u>Q15814</u>
Cytogenetics:	6p21.1
MW:	39.2 kDa