

Product datasheet for **RC203569L3V**

HNRPLL (HNRNPLL) (NM_138394) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HNRPLL (HNRNPLL) (NM_138394) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HNRPLL
Synonyms:	HNRPLL; SRRF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_138394
ORF Size:	1626 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203569).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138394.2
RefSeq Size:	3053 bp
RefSeq ORF:	1629 bp
Locus ID:	92906
UniProt ID:	Q8WVV9
Cytogenetics:	2p22.1
Domains:	RRM
MW:	60.1 kDa


[View online »](#)

Gene Summary:

HNRNPLL is a master regulator of activation-induced alternative splicing in T cells. In particular, it alters splicing of CD45 (PTPRC; MIM 151460), a tyrosine phosphatase essential for T-cell development and activation (Oberdoerffer et al., 2008 [PubMed 18669861]). [supplied by OMIM, Aug 2008]