

Product datasheet for **RC203581**

C11orf49 (NM_001003676) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C11orf49 (NM_001003676) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	C11orf49
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC203581 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGAGTCCGGAGCGCCTAGCCCTACCGACTACGAGTATCTGGCTCAGCGACATGTCCTCACCTACA
TGGAGGATGCAGTGTGCCAGCTGCTAGAAAACAGGAAGATATTAGCCAATATGGAATTGCCAGGTTCTT
CACTGAATATTTAACAGTGTATGCCAGGAACACACATTCTCTTTCGAGAATTCAGCTTCGTCCAAGCC
ACCCCCACAATAGGGTATCATTTCACGGCCCTCTGGAGATGCTTCCGAACTGTGGGCAAAAATGGCG
ATTTGCTGACCATGAAAGAATATCACTGTTTGTGCAATTACTGTGCTCCTGATTTCCCGCTGGAGCTCAC
TCAGAAAAGCAGCCAGGATTGTGCTCATGGACGATGCCATGGACTGCTTGATGCTTTTTTCAGATTTCTC
TTTGCCTTCCAGATCCAGTTTTACTACTCAGAATTCCTGGACAGTGTGGCTGCCATCTATGAGGACCTGC
TGTGAGGCAAGAACCCCAACACAGTATTGTGCCGACGTCGTCCAGTGGGACGACCCGCAACGACCTGC
CTTGGGCGGGCCGGCAGCTGGAGGGCGTGGAGGCGTGGAGGCGTGTCTACCAAGTGTCTGAAAACTGTGT
GATCGGCACAAGTACAGCTGCCACCCAGCACTTGTCAAAGAGGCCCTCAGCAATGTTTCAGAGACTGA
CCTTCTATGGATTCCTCATGGCTCTCTCAAAGCACCGTGAATCAACCAAGCCCTCGGAGCTTTGCCTGA
CAAGGGGATCTGATGCACGACCCAGCAATGGATGAAGAGCTGGAACGGCTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC203581 protein sequence
 Red=Cloning site Green=Tags(s)

MLSPERLALPDYEYLAQRHVLT YMEDAVCQLENREDISQYGIARFFTEYFNSVCQGTHILFREFSFVQA
 TPHNRVSFLRAFWRFCRTV GKNGLL TMKEYHCLLQLLCPDFPLELTQKAARIVLMDDAMDCLMSFSDFL
 FAFQIQFYSEFLDSVAIYEDLLSGKNPNTVIVPTSSSGQHRQRPALGGAGTLEGVEASLFYQCLENLC
 DRHKYSCPPPALVKEALSNVQRLTFYGF LMA LSKHRGINQALGALPDKGDL MHPAMDEELERL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6428_g05.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001003676

ORF Size: 822 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001003676.2](#), [NP_001003676.1](#)

RefSeq Size: 2012 bp

RefSeq ORF: 825 bp


Locus ID: 79096

UniProt ID: [Q9H6J7](#)


Cytogenetics: 11p11.2

MW: 31.3 kDa

Product images:



Circular map for RC203581



Western blot validation of overexpression lysate (Cat# [LY424096]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC203581 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).