

Product datasheet for **RC203584A1V**

Human GMPR2 (NM_016576) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human GMPR2 (NM_016576) AAV Particle
Tag:	Myc-DDK
Symbol:	GMPR2
Synonyms:	GMPR 2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203584 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGACTTCCTGCCTTCCAGCCCTCAGATTCATCGCTACCCGAGGCTAAGCGCCATGCCTCATATTGACA
ACGATGTGAACTGGACTTCAAGGATGTCCTTTTGAGGCCAAACGCAGTACCCTTAAGTCTCGAAGTGA
GGTGGATCTCACAAGATCCTTTTCATTTGGAAGTCAAGCAGACATACTCTGGGGTCCCATCATTGCT
GCCAATATGGATACTGTGGGCACCTTTGAGATGGCCAAGGTTCTCTGTAAGTTCTCTCTTCACTGCTG
TCCATAAGCACTATAGCCTCGTTCAGTGGCAAGAGTTTGCTGGCCAGAATCCTGACTGTCTTGAGCATCT
GGCTGCCAGCTCAGGCACAGGCTCTTCTGACTTTGAGCAGCTGGAACAGATCCTGGAAGCTATTCCCCAG
GTGAAGTATATATGCCTGGATGTGGCAAATGGCTACTCTGAACACTTTGTTGAATTTGTAAGATGTAC
GGAAGCGCTTCCCCAGCACACCATCATGGCAGGGAATGTGGTAACAGGAGAGATGGTAGAAGAGCTCAT
CCTTTCTGGGGCTGACATCATCAAAGTGGGAATTGGGCCAGGCTCTGTGTGTAAGTCTCGAAGAAAAC
GGAGTGGGGTATCCACAGCTCAGCGCAGTGATGGAGTGTGAGATGCTGCTCATGGCCTCAAAGGCCACA
TCATTTAGATGGAGTTGCAGCTGTCTGGGGATGTGGCCAAGGCTTTTGGGGCAGGAGCTGACTTCGT
GATGCTGGGTGGCATGCTGGCTGGGCACAGTGAGTCAGGTGGTGGAGCTCATCGAGAGGGATGGCAAGAAG
TACAAGCTCTTCTATGGAATGAGTTCTGAAATGGCCATGAAGAAGTATGCTGGGGCGTGGCTGAGTACA
GAGCCTCAGAGGAAAGACAGTGAAGTTCTTTTAAAGGAGATGTGGAACATACCATCCGAGACATCCT
AGGAGGGATCCGCTCTACGTGTACCTATGTGGGAGCAGTAAGCTCAAAGAGTTGAGCAGGAGAAGTACC
TTCATCCGAGTCACCCAGCAGGTGAATCCAATCTTCAGTGAGGCGTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC203584 protein sequence Red=Cloning site Green=Tags(s) MTSCLPALRFIATPRLSAMPIDNDVKLDFKDVLLRPKRSTLKSSEVDLTRSFSFRNSKQTYSGVPIIA ANMDTVGTFEMAKVLCKFSLFTAVHKHYSLVQWQEFAGQNPDCLEHLAASSGTGSSDFEQLEQILEAIPQ VKYICLDVANGYSEHFVEFVKDVRKRFPQHTIMAGNVVTGEMVEELILSGADIIVKVGIGPGSVCTTRKKT GVGYPQLSAVMCADAHGLKGHIISDGGCSCPGDVAKAFGAGADFVMLGGMLAGHSESGGELIERDGKK YKLFYGMSSSEMAMKYYAGGVAEYRASEGKTVEVPFKGDVEHTIRDILGGIRSTCTYVGAARKELSRRTT FIRVTQQVNPIFSEAC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_016576
ORF Size:	1098 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_016576.3
RefSeq Size:	1910 bp
RefSeq ORF:	1101 bp
Locus ID:	51292
UniProt ID:	Q9P2T1
Cytogenetics:	14q12
MW:	39.8 kDa