

Product datasheet for **RC203587L4V**

Tryptophanyl tRNA synthetase (WARS) (NM_004184) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tryptophanyl tRNA synthetase (WARS) (NM_004184) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Tryptophanyl tRNA synthetase
Synonyms:	GAMMA-2; HMN9; IFI53; IFP53; WARS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004184
ORF Size:	1413 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203587).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004184.3
RefSeq Size:	2884 bp
RefSeq ORF:	1416 bp
Locus ID:	7453
UniProt ID:	P23381
Cytogenetics:	14q32.2
Domains:	WHEP-TRS, tRNA-synt_1b



[View online »](#)

Protein Families:	Druggable Genome
Protein Pathways:	Aminoacyl-tRNA biosynthesis, Tryptophan metabolism
MW:	53.2 kDa
Gene Summary:	<p>Aminoacyl-tRNA synthetases catalyze the aminoacylation of tRNA by their cognate amino acid. Because of their central role in linking amino acids with nucleotide triplets contained in tRNAs, aminoacyl-tRNA synthetases are thought to be among the first proteins that appeared in evolution. Two forms of tryptophanyl-tRNA synthetase exist, a cytoplasmic form, named WARS, and a mitochondrial form, named WARS2. Tryptophanyl-tRNA synthetase (WARS) catalyzes the aminoacylation of tRNA(trp) with tryptophan and is induced by interferon. Tryptophanyl-tRNA synthetase belongs to the class I tRNA synthetase family. Four transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p>