

Product datasheet for **RC203589A1V**

Human SMAGP (NM_001033873) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human SMAGP (NM_001033873) AAV Particle
Tag: Myc-DDK
Symbol: SMAGP
Synonyms: hSMAGP
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC203589 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGACCAGCCTCCTGACTACTCCTTCTCCAAGAGAAGAACTGATGACCACCCCAATTTTACAGCCCACTG
AGGCCCTGTCCCAGAAGATGGAGCCAGCACAGCACTCATTGCAGTTGTTATCACCGTTGTCTTCCTTAC
CCTGCTCTCGGTGATCTTGTATCTTCTTTACCTGTACAAGAACAAGGCAGCTACGTCACCTATGAA
CCTACAGAAGGTGAGCCAGTGCCATCGTCCAGATGGAGAGTGACTTGGCCAAGGGCAGCGAGAAAGAGG
AATATTTTCATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC203589 protein sequence
 Red=Cloning site Green=Tags(s)

MTSLLTTPSPREELMTTPILQPTTEALSPEDGASTALIAVVITVVFLTLLSVVILIFFYLKYNKGSYVTYE
PTEGEPsAIVQMESDLAKGSEKEYFI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001033873
ORF Size: 291 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001033873.1, NP_001029045.1</u>
RefSeq Size:	1076 bp
RefSeq ORF:	294 bp
Locus ID:	57228
UniProt ID:	<u>Q0VAQ4</u>
Cytogenetics:	12q13.13
MW:	10.7 kDa