

Product datasheet for **RC203612L4V**

ERGIC1 (NM_001031711) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ERGIC1 (NM_001031711) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ERGIC1
Synonyms:	AMC2; AMCN; ERGIC-32; ERGIC32; NET24
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001031711
ORF Size:	870 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203612).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001031711.1
RefSeq Size:	2949 bp
RefSeq ORF:	873 bp
Locus ID:	57222
UniProt ID:	Q969X5
Cytogenetics:	5q35.1
Protein Families:	Transmembrane
MW:	32.6 kDa



[View online »](#)

Gene Summary:

This gene encodes a cycling membrane protein which is an endoplasmic reticulum-golgi intermediate compartment (ERGIC) protein which interacts with other members of this protein family to increase their turnover. [provided by RefSeq, Jul 2008]