

## Product datasheet for **RC203613L4V**

### **U2AF65 (U2AF2) (NM\_001012478) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	U2AF65 (U2AF2) (NM_001012478) Human Tagged ORF Clone Lentiviral Particle
Symbol:	U2AF2
Synonyms:	U2AF65
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001012478
ORF Size:	1413 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203613).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001012478.1</a> , <a href="#">NP_001012496.1</a>
RefSeq Size:	3136 bp
RefSeq ORF:	1416 bp
Locus ID:	11338
UniProt ID:	<a href="#">P26368</a>
Cytogenetics:	19q13.42
Protein Pathways:	Spliceosome
MW:	52.9 kDa



[View online »](#)

**Gene Summary:**

U2 auxiliary factor (U2AF), comprised of a large and a small subunit, is a non-snRNP protein required for the binding of U2 snRNP to the pre-mRNA branch site. This gene encodes the U2AF large subunit which contains a sequence-specific RNA-binding region with 3 RNA recognition motifs and an Arg/Ser-rich domain necessary for splicing. The large subunit binds to the polypyrimidine tract of introns early during spliceosome assembly. Multiple transcript variants have been detected for this gene, but the full-length nature of only two have been determined to date. [provided by RefSeq, Jul 2008]