

Product datasheet for **RC203637L4V**

RDH13 (NM_138412) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RDH13 (NM_138412) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RDH13
Synonyms:	SDR7C3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138412
ORF Size:	780 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203637).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138412.2
RefSeq Size:	2023 bp
RefSeq ORF:	783 bp
Locus ID:	112724
UniProt ID:	Q8NBN7
Cytogenetics:	19q13.42
Domains:	adh_short
Protein Families:	Druggable Genome



[View online »](#)

MW: 28.8 kDa

Gene Summary: This gene encodes a mitochondrial short-chain dehydrogenase/reductase, which catalyzes the reduction and oxidation of retinoids. The encoded enzyme may function in retinoic acid production and may also protect the mitochondria against oxidative stress. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2009]