

Product datasheet for **RC203647L3V**

SETDB1 (NM_012432) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SETDB1 (NM_012432) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SETDB1
Synonyms:	ESET; H3-K9-HMTase4; KG1T; KMT1E; TDRD21
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012432
ORF Size:	3870 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203647).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012432.2 , NP_036564.2
RefSeq Size:	4446 bp
RefSeq ORF:	3873 bp
Locus ID:	9869
UniProt ID:	Q15047
Cytogenetics:	1q21.3
Domains:	SET, MBD, TUDOR, PreSET, Pre-SET
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Lysine degradation

MW: 143 kDa

Gene Summary: This gene encodes a histone methyltransferase which regulates histone methylation, gene silencing, and transcriptional repression. This gene has been identified as a target for treatment in Huntington Disease, given that gene silencing and transcription dysfunction likely play a role in the disease pathogenesis. Alternatively spliced transcript variants of this gene have been described.[provided by RefSeq, Jun 2011]