

## Product datasheet for **RC203648L2V**

### TM9SF4 (NM\_014742) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TM9SF4 (NM_014742) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TM9SF4
Synonyms:	dj836N17.2
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_014742
ORF Size:	1875 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203648).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_014742.3</a> , <a href="#">NP_055557.1</a>
RefSeq Size:	3996 bp
RefSeq ORF:	1929 bp
Locus ID:	9777
UniProt ID:	<a href="#">Q92544</a>
Cytogenetics:	20q11.21
Domains:	EMP70
Protein Families:	Transmembrane



[View online »](#)

**MW:** 72.54 kDa

**Gene Summary:** Associates with proteins harboring glycine-rich transmembrane domains and ensures their efficient localization to the cell surface (PubMed:25999474). Regulates the assembly and activity of V-ATPase in colon cancer cells via its interaction with V-type proton ATPase subunit H (ATP6V1H) and contributes to V-ATPase-mediated pH alterations in cancer cells which play an important role in drug resistance and invasiveness of colon cancer cells (PubMed:25659576). Plays an important role in an atypical phagocytic activity of metastatic melanoma cells called cannibalism and is involved in the pH regulation of the intracellular vesicles in tumor cells (PubMed:19893578).[UniProtKB/Swiss-Prot Function]