

Product datasheet for **RC203664L3V**

C21orf70 (FAM207A) (NM_058190) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	C21orf70
Synonyms:	C21orf70; PRED56
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_058190
ORF Size:	690 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC203664).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_058190.2
RefSeq Size:	964 bp
RefSeq ORF:	693 bp
Locus ID:	85395
UniProt ID:	Q9NSI2
Cytogenetics:	21q22.3



MW: 25.5 kDa