

Product datasheet for **RC203670L4V**

C14orf80 (NM_173608) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C14orf80 (NM_173608) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C14orf80
Synonyms:	chromosome 14 open reading frame 80; MGC16771
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_173608
ORF Size:	744 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203670).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_173608.1 , NP_775879.1
RefSeq Size:	1563 bp
RefSeq ORF:	1124 bp
Locus ID:	283643
Cytogenetics:	14q32.33
MW:	27.7 kDa
Gene Summary:	May play a role in counteracting perturbation of actin filaments, such as after treatment with the actin depolymerizing microbial metabolite Chivosazole F.[UniProtKB/Swiss-Prot Function]



[View online »](#)