

Product datasheet for **RC203677L3V**

SLC35F5 (NM_025181) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Symbol: SLC35F5

Mammalian Cell Selection: Puromycin

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Tag: Myc-DDK

ACCN: NM_025181

ORF Size: 1569 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC203677).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_025181.2](#)

RefSeq Size: 4412 bp

RefSeq ORF: 1572 bp

Locus ID: 80255

UniProt ID: [Q8WV83](#)

Cytogenetics: 2q14.1

Domains: DUF6



Protein Families: Transmembrane

MW: 58.9 kDa

Gene Summary: Putative solute transporter.[UniProtKB/Swiss-Prot Function]