

Product datasheet for **RC203679**

EPS8 like protein 3 (EPS8L3) (NM_133181) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EPS8 like protein 3 (EPS8L3) (NM_133181) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	EPS8 like protein 3
Synonyms:	EPS8R3; HYPT5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC203679 representing NM_133181.
 Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAAATTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTCAAGGCCAGCAGCAGAGCCATTTACTTGCACCGGAAGGAGTACTCCAGAACCTCACCTCAGAG
CCCACCTCCTGCAGCACAGGGTGGAGCACTTGATGACATGCAAGCAGGGGAGTCAGAGAGTCCAGGGG
CCCGAGGATGCCTTGCAGAAAGCTGTTTCGAGATGGATGCACAGGGCCGGGTGTGGAGCCAAGACTTGATC
CTGCAGGTGAGGACGGCTGGCTGCAGTCTGCTGGACATTGAGACCAAGGAGGAGTGGACTCTTACCGC
CTAGACAGCATCCAGGCCATGAATGTGGCGCTCAACACATGTTCTACAACCTCCATCCTGTCCATCACC
GTGCAGGAGCCGGGCTGCCAGGCACTAGCACTCTGCTCTTCCAGTGCCAGGAAGTGGGGGAGAGCGA
CTGAAGACCAGCTGCAGAAAGCTCTGGAGGAAGAGCTGGAGCAAAGACCTCGACTTGGAGGCTTCAG
CCAGGCCAGGACAGATGGAGGGGCTGCTATGGAAAGGCCGCTCCCTATGGAGCAGGCACGCTATCTG
GAGCCGGGATCCCTCCAGAACAGCCCCACCAGAGGACCTAGAGCACAGCTCCACCACCCCCAAGG
CCCCTGCCAGCCACACCACTGCCCCGAGAACCAAGTGCCTTTACTCTGCCTCCTCAAGGCGGTCTCT
TCCCCCAGGACCCAGAGAGGGACGAGGAAGTGTGAACCATGTCTAAGGGACATTGAGCTGTTATG
GGAAAGCTGGAGAAGGCCAGGCAAGACAGCAGGAAGAAGAAATTTGGAAAAAACAAGGACCAG
GGAGGTCTACCCAGGCACAGTACATTGACTGCTTCCAGAAGTCAAGCACAGCTTCAACCTCCTGGGA
AGGCTGGCCACCTGGCTGAAGGAGACAAGTGCCCTGAGCTCGTACACATCCTCTTCAAGTCCCTGAAC
TTCATCCTGGCCAGGTGCCCTGAGGCTGGCTAGCAGCCCAAGTATCTACCCCTCCTACCCCTAAA
GCTATCAACCTGCTACAGTCTGTCTAAGCCACCTGAGAGTAACCTTTGGATGGGGTGGGGCCAGCC
TGGACCACTAGCCGGGCGACTGGACAGCGATGAGCCCTGCCCTACCAACCCACATTCAGATGAC
TGGCAACTTCCAGAGCCCTCCAGCCAAGCACCTTAGGATACCAGGACCTGTTTCCCTTCGGCGGGGA
AGTCATAGGTAGGGAGCACCTCACATTTCTCAGGAGAAGACACACAACCATGACCTCAGCCTGGG
GACCCCAACTCCAGGCCCTCCAGCCCCAACCTGCCAGCCAGCCCTGAAAATGCAAGTCTTGACGAG
TTTGAAGCTAGGAACCCACGGGAAGTACTGTGGTCCAGGGAGAGAAGCTGGAGGTTCTGGACCACAGC
AAGCGGTGGTGGTGGTGAAGAATGAGGCGGACGGAGCGGTACATTCCAAGCAACATCTGGAGCCC
CTACAGCCGGGACCCCTGGGACCCAGGGCCAGTACCCTCTCGGGTTCCATGCTTCGACTTAGCTCG
AGGCTGAAGAGGTCACAGACTGGCTGCAGGCAGAGAATTCTCCACTGCCACGGTGAGGACACTGGG
TCCCTGACGGGAGCCAGCTACTTCGCATAAGACCTGGGAGCTACAGATGCTATGTCCACAGGAGGCC
CCACGAATCCTGTCCCGCTGGAGGCTGTCAGAAGGATGCTGGGGATAAGCCCT
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```


Protein Sequence:

>Peptide sequence encoded by RC203679
 Blue=ORF Red=Cloning site Green=Tag(s)

```
MSRPSSRAIYLHRKEYSQNLTSPTLLQHRVEHLMTCQKQSQRVQGPEDALQKLFEQDAQGRVWSQDLI
LQVRDGLQLLDIETKEELDSYRLDSIQAMNVALNTCSYNSILSITVQEPGLPGTSTLLFQCEVGAER
LKTSLQKALEEELEQRPRLLGGLQPGQDRWRGPAMERPLPMEQARYLEPGIPPEQPHQRTLEHSLPPSPR
PLPRHTSAREPSAFTLPPPRSSPEDPERDEEVLNHVLRDIELFMGKLEKAQAKTSRKKKFGKKNKDQ
GGLTQAQYIDCFQKIKHSFNLLGRLATWLKETSAPELVHILFKSLNFILARCPEAGLAAQVISPLTPK
AINLLQSCLSPPESNLWMGLGPAWTTSRADWTGDEPLPYQPTFSDDWQLPEPSSQAPLGYQDPVSLRRG
SHRLGSTSHFPQEKTHNDPQPGDPNSRPSPKPAQPALKMQVL YEFARNPRELTVVQGEKLEVL DHS
KRWWLVKNEAGRSGYIPSNILEPLQPGTGTGQSPSRVPMRLRSSRPEEVDWLQAENFSTATVRTLG
SLTGSQLLRIRPGLQMLCPQEAPRILSRLEAVRRMLGISP
TRTRPLEQKLISEEDLAANDILDYKDDDDKV
```

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_133181

ORF Size: 1779 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_133181.4](#)


RefSeq Size: 2329 bp

RefSeq ORF: 1782 bp

Locus ID: 79574

UniProt ID:	<u>Q8TE67</u>
Cytogenetics:	1p13.3
Protein Families:	Druggable Genome
MW:	66.9 kDa
Gene Summary:	This gene encodes a protein that is related to epidermal growth factor receptor pathway substrate 8 (EPS8), a substrate for the epidermal growth factor receptor. The function of this protein is unknown. Alternatively spliced transcript variants encoding different isoforms exist. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC203679