

Product datasheet for **RC203679L3V**

EPS8 like protein 3 (EPS8L3) (NM_133181) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	EPS8 like protein 3 (EPS8L3) (NM_133181) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EPS8 like protein 3
Synonyms:	EPS8R3; HYPT5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_133181
ORF Size:	1779 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203679).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133181.3
RefSeq Size:	2329 bp
RefSeq ORF:	1782 bp
Locus ID:	79574
UniProt ID:	Q8TE67
Cytogenetics:	1p13.3
Protein Families:	Druggable Genome
MW:	66.9 kDa


[View online »](#)

Gene Summary:

This gene encodes a protein that is related to epidermal growth factor receptor pathway substrate 8 (EPS8), a substrate for the epidermal growth factor receptor. The function of this protein is unknown. Alternatively spliced transcript variants encoding different isoforms exist. [provided by RefSeq, Jul 2008]