

Product datasheet for **RC203685L4V**

ECHDC3 (NM_024693) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ECHDC3 (NM_024693) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ECHDC3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024693
ORF Size:	909 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203685).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024693.2
RefSeq Size:	1674 bp
RefSeq ORF:	912 bp
Locus ID:	79746
UniProt ID:	Q96DC8
Cytogenetics:	10p14
Domains:	ECH
MW:	32.7 kDa



[View online »](#)

Gene Summary: May play a role in fatty acid biosynthesis and insulin sensitivity.[UniProtKB/Swiss-Prot Function]