

## Product datasheet for **RC203689L4V**

### UBE2Z (NM\_023079) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	UBE2Z (NM_023079) Human Tagged ORF Clone Lentiviral Particle
Symbol:	UBE2Z
Synonyms:	HOYS7; USE1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_023079
ORF Size:	738 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203689).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_023079.3</a> , <a href="#">NP_075567.1</a>
RefSeq Size:	3101 bp
RefSeq ORF:	1065 bp
Locus ID:	65264
UniProt ID:	<a href="#">Q9H832</a>
Cytogenetics:	17q21.32
Domains:	UBCc
Protein Pathways:	Acute myeloid leukemia, Pathways in cancer, Ubiquitin mediated proteolysis



[View online »](#)

**MW:** 38 kDa

**Gene Summary:** This gene encodes an enzyme which ubiquitinates proteins which participate in signaling pathways and apoptosis. [provided by RefSeq, Feb 2012]