

Product datasheet for **RC203706L3V**

GBA2 (NM_020944) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GBA2 (NM_020944) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GBA2
Synonyms:	AD035; NLGase; SPG46
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_020944
ORF Size:	2781 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203706).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020944.2
RefSeq Size:	3639 bp
RefSeq ORF:	2784 bp
Locus ID:	57704
UniProt ID:	Q9HCG7
Cytogenetics:	9p13.3
Domains:	DUF608
Protein Families:	Druggable Genome



[View online »](#)

MW: 104.6 kDa

Gene Summary: This gene encodes a microsomal beta-glucosidase that catalyzes the hydrolysis of bile acid 3-O-glucosides as endogenous compounds. Studies to determine subcellular localization of this protein in the liver indicated that the enzyme was mainly enriched in the microsomal fraction where it appeared to be confined to the endoplasmic reticulum. This putative transmembrane protein is thought to play a role in carbohydrate transport and metabolism. [provided by RefSeq, Jul 2008]