

Product datasheet for **RC203722L4V**

OTUD6B (NM_016023) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	OTUD6B (NM_016023) Human Tagged ORF Clone Lentiviral Particle
Symbol:	OTUD6B
Synonyms:	CGI-77; DUBA-5; DUBA5; IDDFSDA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016023
ORF Size:	879 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203722).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016023.3 , NP_057107.2
RefSeq Size:	3306 bp
RefSeq ORF:	882 bp
Locus ID:	51633
UniProt ID:	Q8N6M0
Cytogenetics:	8q21.3
Domains:	OTU
Protein Families:	Protease



[View online »](#)

MW: 33.8 kDa

Gene Summary: This gene encodes a member of the ovarian tumor domain (OTU)-containing subfamily of deubiquitinating enzymes. Deubiquitinating enzymes are primarily involved in removing ubiquitin from proteins targeted for degradation. This protein may function as a negative regulator of the cell cycle in B cells. [provided by RefSeq, Nov 2013]