

Product datasheet for **RC203737A1V**

Human UCHL5IP (HAUS7) (NM_017518) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human UCHL5IP (HAUS7) (NM_017518) AAV Particle
Tag:	Myc-DDK
Symbol:	UCHL5IP
Synonyms:	UCHL5IP; UIP1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203737 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGCGGGGCGGGCTCGGAGCGCGAAACATGGCGGGGCGAGACGCTGGCTGCGGCCGTGGCGGCGACG
ACTACTCAGAGGACGAGGGCGACAGCAGCGTGTCCAGGGCGGCTGTGGAGGTGTTCCGGGAAGCTGAAGGA
CCTAAACTGCCCTTCTCGAGGGTCTGTATATCACAGAGCCAAAGACAATTAGGAAGTGTGTGCAGC
CCCTCAGAGTACCGCTTGGAGATCCTAGAGTGGATGTGTACCCGGGTCTGGCCCTCACTGCAGGACAGGT
TCAGTCACTGAAAGGGGTCCCAACAGAGGTGAAGATCCAAGAAATGACGAAGCTGGGCCACGAGCTGAT
GCTGTGTGCGCCAGATGACCAGGAGTCCCTCAAGGGCTGTGCCTGCGCCCAAGCAGCTACACTTCATG
GACCAAGTTGCTCGATACCATCCGAGCCTGACCATTGGGTGCTCCAGTTGCTCGAGCCTGATGGAGCACT
TCGAGGACACCAGGGAGAAGAACGAGGCCTTGCTGGGGGAGCTCTTCTAGCCCCACCTGCAGATGCT
CCTGAATCCAGAGTGCAGCCGTGGCCCTGGACATGCAGCCCTCCTCAACAAGCAGAGTGATGACTGG
CAGTGGGCCAGTGCCTCTGCCAAGTCCGAGGAGGAGAGAAGCTGGCGGAGCTTGCCAGGCAGCTGCAGG
AGAGTGCTGCCAAGTTGCACGCGCTTAGAACGGAGTACTTTGCACAGCATGAGCAAGGGGCTGCTGCGGG
CGCAGCCGACATCAGCACCTAGACCAGAAGCTGCGTCTGGTCACTTCCGACTTCCACCAGCTAATCTTG
GCTTTTCTCCAAGTCTACGACGACGAGCTGGGCGAGTGCTGCCAGCGCCAGGCCCTGACCTCCACCCGT
GCGGCCCATCATCCAGGCCACGCACCAGAATCTGACTTCCTACAGCCAACTGCTGCAAGTGGTCATGGC
AGTTGCTGACACCTCTGCGAAGGCCGTGGAGACCGTGAAGAAGCAGCAAGGCGAGCAGATCTGCTGGGGT
GGCAGCAGCTCCGTATGAGTCTAGCTACCAAGATGAATGAACCTAATGGAGAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC203737 protein sequence Red=Cloning site Green=Tags(s) MGGARLGARNMAGQDAGCGRGDDYSEDEGDSSVSRAAVEVFGKLDLNCPLFLEGLYITEPKTIQELLCS PSEYRLEILEWMCTRVWPSLQDRFSSLKGVPTVEVKIQEMTKLGHEMLCAPDDQELLKGCACAQKQLHFM DQLLDTIRSLTIGCSSCSLMEHFEDTREKNEALLGELFSSPHLQMLLNPECDPWPLDMQPLLNKQSDDW QWASASAKSEEEEKLAELARQLQESA AKLHALRTEYFAQHEQGAAAGAADISTLDQKLRLVTSDFHQLIL AFLQVYDDELGECCQRPGLHPCGPIIQATHQNLTSYSQLLQVMAVADTSAKAVETVKKQQGEQICWG GSSSVMSLATKMNELMEK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_017518
ORF Size:	1104 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_017518.5
RefSeq Size:	1378 bp
RefSeq ORF:	1107 bp
Locus ID:	55559
UniProt ID:	Q99871
Cytogenetics:	Xq28
MW:	40.8 kDa