

Product datasheet for **RC203737L4V**

UCL5IP (HAUS7) (NM_017518) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	UCL5IP (HAUS7) (NM_017518) Human Tagged ORF Clone Lentiviral Particle
Symbol:	UCL5IP
Synonyms:	UCL5IP; UIP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017518
ORF Size:	1104 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203737).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017518.5
RefSeq Size:	1378 bp
RefSeq ORF:	1107 bp
Locus ID:	55559
UniProt ID:	Q99871
Cytogenetics:	Xq28
MW:	40.8 kDa


[View online »](#)

Gene Summary:

This gene encodes a subunit of the augmin complex, which regulates centrosome and mitotic spindle integrity, and is necessary for the completion of cytokinesis. The encoded protein was identified by interaction with ubiquitin C-terminal hydrolase 37. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2012]