

## Product datasheet for **RC203739L4V**

### TMEM59L (NM\_012109) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TMEM59L (NM_012109) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMEM59L
Synonyms:	BSMAP; C19orf4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012109
ORF Size:	1026 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203739).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012109.2</a>
RefSeq Size:	1712 bp
RefSeq ORF:	1029 bp
Locus ID:	25789
UniProt ID:	<a href="#">Q9UK28</a>
Cytogenetics:	19p13.11
Protein Families:	Transmembrane
MW:	37.6 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a predicted type-I membrane glycoprotein. The encoded protein may play a role in functioning of the central nervous system. [provided by RefSeq, Jul 2008]