

Product datasheet for **RC203746A1V**

Human U1SNRNPBP (SNRNP35) (NM_180699) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human U1SNRNPBP (SNRNP35) (NM_180699) AAV Particle
Tag:	Myc-DDK
Symbol:	U1SNRNPBP
Synonyms:	HM-1; U1SNRNPBP
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203746 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAAACCAGCTAACATGAACGATTGGATGCCCATCGCCAAGGAGTATGATCCACTCAAAGCGGGCAGCA
TTGATGGCACCAGATGAAGACCCACACGACCGCGCGGTCTGGAGGGCAATGCTGGCAGCATATGTCCCAA
CAAAGGTGTCATAGGAGATCCCCTCCTCACCCTGTTTGTGGCCAGACTAACTTGCAGACCAAGGAGGAC
AAATTAAGGAAGTCTTTCCCGCTATGGTGACATCCGGCGGCTTCGGCTGGTCAGGGACTTGGTCACAG
GTTTTTCAAAGGGCTACGCCTTCATCGAATACAAGGAGGAGCGTGCCGTGATCAAAGCTTACCGAGATGC
TGATGGCCTGGTTATTGACCAGCATGAGATATTTGTGGACTACGAGCTGGAAAGGACTCTCAAAGGGTGG
ATCCCTCGGCGACTTGGAGGCGGTCTTGGGGGAAAAAAGGAGTCTGGGCAACTGAGATTTGGGGGACGGG
ACCGGCCTTTTCGAAAACCTATTAACCTTGCCAGTTGTTAAAAACGACCTCTATAGAGAGGGAAAAACGGGA
AAGGCGGGAGCGATCTCGATCCCGAGAAAGACACTGGGACTCGAGGACAAGGGATCGAGACCATGACAGG
GGCCGGGAGAAGAGATGGCAAGAAAGAGAGCCGACCAGGGTGTGGCCCGACAATGACTGGGAGAGAGAGA
GGGACTTCAGAGATGACAGGATCAAGGGGAGGGAGAAGAAGGAAAGAGGCAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence:	>RC203746 protein sequence Red=Cloning site Green=Tags(s) MKPANMNDWMPiAKEYDPLKAGSIDGTDEDPHDRAVWRAMLARYVPNKGVIQDPLLTLFVARLNLQTKED KLKEVFSRYGDIRRLRLVRDLVTGFSKGYAFIEYKEERAVIKAYRDADGLVIDQHEIFVDYELERTLKGW IPRRLGGGLGGKKEGQLRFGRDRPFRKPINLPVVKNDLYREGKRERRERSRSRERHWDSTRDRDHDR GREKRWQEREPTRVWPDNDWERERDFRDDRIGREKKERK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_180699
ORF Size:	753 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_180699.1, NP_851030.1</u>
RefSeq Size:	1070 bp
RefSeq ORF:	756 bp
Locus ID:	11066
UniProt ID:	<u>Q16560</u>
Cytogenetics:	12q24.31
MW:	30 kDa