

Product datasheet for **RC203751L1V**

VSIG4 (NM_007268) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | VSIG4 (NM_007268) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | VSIG4 |
| Synonyms: | CRlg; Z39IG |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_007268 |
| ORF Size: | 1197 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC203751). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_007268.1 |
| RefSeq Size: | 1869 bp |
| RefSeq ORF: | 1200 bp |
| Locus ID: | 11326 |
| UniProt ID: | Q9Y279 |
| Cytogenetics: | Xq12 |
| Domains: | ig, IGv, IGc2, IG |
| Protein Families: | Transmembrane |



[View online »](#)

MW: 44 kDa

Gene Summary: This gene encodes a v-set and immunoglobulin-domain containing protein that is structurally related to the B7 family of immune regulatory proteins. The encoded protein may be a negative regulator of T-cell responses. This protein is also a receptor for the complement component 3 fragments C3b and iC3b. Alternate splicing results in multiple transcript variants. [provided by RefSeq, May 2010]