

Product datasheet for **RC203764L4V**

OXCT1 (NM_000436) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	OXCT1 (NM_000436) Human Tagged ORF Clone Lentiviral Particle
Symbol:	OXCT1
Synonyms:	OXCT; SCOT
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_000436
ORF Size:	1560 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203764).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000436.3
RefSeq Size:	3572 bp
RefSeq ORF:	1563 bp
Locus ID:	5019
UniProt ID:	P55809
Cytogenetics:	5p13.1
Domains:	CoA_trans



[View online »](#)

Protein Pathways:	Butanoate metabolism, Synthesis and degradation of ketone bodies, Valine, leucine and isoleucine degradation
MW:	56.2 kDa
Gene Summary:	This gene encodes a member of the 3-oxoacid CoA-transferase gene family. The encoded protein is a homodimeric mitochondrial matrix enzyme that plays a central role in extrahepatic ketone body catabolism by catalyzing the reversible transfer of coenzyme A from succinyl-CoA to acetoacetate. Mutations in this gene are associated with succinyl CoA:3-oxoacid CoA transferase deficiency. [provided by RefSeq, Jul 2008]