

Product datasheet for RC203767L3

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US
Phone: +1-888-267-4436
https://www.origene.com
techsupport@origene.com
EU: info-de@origene.com
CN: techsupport@origene.cn

SERPING1 (NM_000062) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: SERPING1 (NM_000062) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: SERPING1

Synonyms: C1IN; C1INH; C1NH; HAE1; HAE2

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC203767).

Sequence:

Sgfl-Mlul

Restriction Sites: Cloning Scheme:





 $[\]ensuremath{^*}$ The last codon before the Stop codon of the ORF.

ACCN: NM_000062

ORF Size: 1500 bp





SERPING1 (NM_000062) Human Tagged Lenti ORF Clone - RC203767L3

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 000062.2</u>

RefSeq Size: 1984 bp RefSeq ORF: 1503 bp Locus ID: 710

UniProt ID: P05155

Cytogenetics: 11q12.1

Domains: SERPIN

Protein Families: Druggable Genome, Secreted Protein

Protein Pathways: Complement and coagulation cascades

MW: 55.2 kDa


Gene Summary: This gene encodes a highly glycosylated plasma protein involved in the regulation of the

complement cascade. Its encoded protein, C1 inhibitor, inhibits activated C1r and C1s of the first complement component and thus regulates complement activation. It is synthesized in the liver, and its deficiency is associated with hereditary angioneurotic oedema (HANE). Alternative splicing results in multiple transcript variants encoding the same isoform.

[provided by RefSeq, May 2020]



Product images:



Circular map for RC203767L3