

Product datasheet for RC203776L4

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CDw75 (ST6GAL1) (NM_173216) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: CDw75 (ST6GAL1) (NM_173216) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: CDw75

Synonyms: SIAT1; ST6Gall; ST6N

Mammalian Cell Puromycin

Selection: Vector:

pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clor

Sequence:

The ORF insert of this clone is exactly the same as(RC203776).

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_173216

ORF Size: 1218 bp





CDw75 (ST6GAL1) (NM_173216) Human Tagged Lenti ORF Clone - RC203776L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 173216.1</u>

 RefSeq Size:
 4604 bp

 RefSeq ORF:
 1221 bp

 Locus ID:
 6480

 UniProt ID:
 P15907

Cytogenetics: 3q27.3

Protein Families: Secreted Protein

Protein Pathways: Metabolic pathways, N-Glycan biosynthesis

MW: 46.6 kDa

Gene Summary: This gene encodes a member of glycosyltransferase family 29. The encoded protein is a type

II membrane protein that catalyzes the transfer of sialic acid from CMP-sialic acid to


galactose-containing substrates. The protein, which is normally found in the Golgi but can be proteolytically processed to a soluble form, is involved in the generation of the cell-surface carbohydrate determinants and differentiation antigens HB-6, CD75, and CD76. This gene has

been incorrectly referred to as CD75. Alternative splicing results in multiple transcript


variants. [provided by RefSeq, Jul 2017]



Product images:



Circular map for RC203776L4



Double digestion of RC203776L4 using Sgfl and Mlul