

Product datasheet for RC203789A1V

Human MCP4 (CCL13) (NM_005408) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human MCP4 (CCL13) (NM_005408) AAV Particle
Tag: Myc-DDK
Symbol: MCP4
Synonyms: CKb10; MCP-4; NCC-1; NCC1; SCYA13; SCYL1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC203789 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGAAAGTCTCTGCAGTGCTTCTGTGCCTGCTGCTCATGACAGCAGCTTTCAACCCCAAGGACTTGCTC
 AGCCAGATGCACTCAACGTCCCATCTACTTGCTGCTTCACATTTAGCAGTAAGAAGATCTCCTTGACAG
 GCTGAAGAGCTATGTGATCACCACCAGCAGGTGTCCCAAGGCTGTCATCTTCAGAACCAAACTGGGC
 AAGGAGATCTGTGCTGACCCAAAGGAGAAGTGGGTCCAGAATTATATGAAACACCTGGGCCGAAAGCTC
 ACACCTGAAGACT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC203789 protein sequence
 Red=Cloning site Green=Tags(s)
 MKVSAVLLCLLLMTAAFPQGLAQPDALNVPSTCCFTFSSKKISLQRLKSYVITTSRCPQKAVIFRTKLG
 KEICADPKEKWVQNYMKHLGRKAHTLKT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_005408
ORF Size: 294 bp


[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_005408.2
RefSeq Size:	861 bp
RefSeq ORF:	297 bp
Locus ID:	6357
UniProt ID:	Q99616
Cytogenetics:	17q12
MW:	11 kDa