

Product datasheet for **RC203822A1V**

Human ZFAND2B (NM_138802) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ZFAND2B (NM_138802) AAV Particle
Tag:	Myc-DDK
Symbol:	ZFAND2B
Synonyms:	AIRAPL
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203822 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGTTTCCGGACCTCGGCGCTCACTGTTCCGAGCCGAGCTGTCAGCGCTTGGATTTCTGCCGCTTA
AGTGTGATGCCTGCTCAGGCATCTTCTGCGCAGACCATGTGGCTACGCCCAGCATCACTGTGGATCTGC
TTACAAAAGGATATCCAGGTACCTGTGTGCCCTCTCTGTAATGTGCCTGTGCCTGTGGCCAGAGGGGAG
CCCCTGACCGTGCTGTGGGAGAGCACATTGACAGAGACTGTCGCTCTGATCCAGCACAGCAAAACGTA
AGATCTTCACCAATAAGTGTGAACGCGCTGGCTGCCGCGAGGAAATGATGAACTGACCTGTGAACG
CTGTAGCCGAACTTCTGCATCAAGCACCGGCATCCACTGGACCATGATTGCTCTGGGGAGGGGCACCCA
ACCAGCCGGGCAGGACTTGCTGCCATCTCCAGAGCACAAAGCTGTGGCTTCTACAAGCACTGTCCCCAGCC
CAAGTCAAACCATGCCTTCTGTACCTCTCCAGCAGAGCCACAACCCGATCTCCGTCTGGACAGCCCC
TCCAGTGATTGCTTGCAGAATGGCCTGAGTGAGGATGAAGCTCTGCAGCGGGCCCTGGAAATGTCCCTG
GCAGAAACCAACCCAGGTTCCAAGTTGTGAGGAGGAAGAAGACCTAGCTTTAGCACAAGCACTGTGAG
CCAGTGAGGCAGAATACCAGCGGCAGCAGGCCAGAGCCGAGCTCGAAGCCGTCCAACCTGCAGCCTGTG
C

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC203822 protein sequence Red=Cloning site Green=Tags(s) MEFPDLGAHCSEPCQRLDFLPLKCDACSGIFCADHVAYAQHHCQSAYQKDIQVPVCPLCNVPVPVARGE PPDRAVGEHIDRDCRSDPAQQKRKIFTNKCERAGCRQREMMKLTCECSRNFCEIKHRHPLDHDCSGEGHP TSRAGLAASRAQAVASTSTVSPSQTMPSTSPSRATTRSPSWTAPPVIALQNGLSEDEALQRALEMSL AETKPQVPSCQEEEDLALAQAALSAEYQRQQAQSRSSKPSNCSLC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_138802
ORF Size:	771 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_138802.1</u>
RefSeq Size:	1323 bp
RefSeq ORF:	774 bp
Locus ID:	130617
UniProt ID:	<u>Q8WV99</u>
Cytogenetics:	2q35
MW:	28 kDa