

## Product datasheet for **RC203822L4V**

### ZFAND2B (NM\_138802) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ZFAND2B (NM_138802) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZFAND2B
Synonyms:	AIRAPL
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138802
ORF Size:	771 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203822).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_138802.1</a>
RefSeq Size:	1323 bp
RefSeq ORF:	774 bp
Locus ID:	130617
UniProt ID:	<a href="#">Q8WV99</a>
Cytogenetics:	2q35
Domains:	ZnF_AN1, UIM
MW:	28 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein containing AN1-type zinc-fingers and ubiquitin-interacting motifs. The encoded protein likely associates with the proteasome to stimulate the degradation of toxic or misfolded proteins. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Aug 2012]