

Product datasheet for **RC203832L3V**

C1QC (NM_172369) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | C1QC (NM_172369) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | C1QC |
| Synonyms: | C1Q-C; C1QG |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_172369 |
| ORF Size: | 735 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC203832). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_172369.2 |
| RefSeq Size: | 1182 bp |
| RefSeq ORF: | 738 bp |
| Locus ID: | 714 |
| UniProt ID: | P02747 |
| Cytogenetics: | 1p36.12 |
| Protein Families: | Secreted Protein |
| Protein Pathways: | Complement and coagulation cascades, Prion diseases, Systemic lupus erythematosus |



[View online »](#)

MW: 25.8 kDa

Gene Summary: This gene encodes the C-chain polypeptide of serum complement subcomponent C1q, which associates with C1r and C1s to yield the first component of the serum complement system. C1q is composed of 18 polypeptide chains which include 6 A-chains, 6 B-chains, and 6 C-chains. Each chain contains an N-terminal collagen-like region and a C-terminal C1q globular domain. C1q deficiency is associated with lupus erythematosus and glomerulonephritis. [provided by RefSeq, Dec 2016]