

Product datasheet for **RC203837L4V**

FAM114A1 (NM_138389) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FAM114A1 (NM_138389) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FAM114A1
Synonyms:	Noxp20
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138389
ORF Size:	1689 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203837).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138389.1 , NP_612398.1
RefSeq Size:	4138 bp
RefSeq ORF:	1692 bp
Locus ID:	92689
UniProt ID:	Q8IWE2
Cytogenetics:	4p14
MW:	60.8 kDa


[View online »](#)

Gene Summary:

The protein encoded by this gene belongs to the FAM114 family and may play a role in neuronal cell development. Alternative splicing results in multiple transcript variants.
[provided by RefSeq, Apr 2017]