

Product datasheet for **RC203842L1V**

Granuphilin (SYTL4) (NM_080737) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Granuphilin (SYTL4) (NM_080737) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Granuphilin
Synonyms:	SLP4
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_080737
ORF Size:	2013 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203842).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_080737.1 , NP_542775.1
RefSeq Size:	3936 bp
RefSeq ORF:	2016 bp
Locus ID:	94121
UniProt ID:	Q96C24
Cytogenetics:	Xq22.1
Domains:	C2, RPH3A_effector
MW:	76 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the synaptotagmin like protein family. Members of this family are characterized by an N-terminal Rab27 binding domain and C-terminal tandem C2 domains. The encoded protein binds specific small Rab GTPases and is involved in intracellular membrane trafficking. This protein binds Rab27 and may be involved in inhibiting dense core vesicle exocytosis. Alternate splicing results in multiple transcript variants that encode the same protein. [provided by RefSeq, Mar 2010]