

Product datasheet for RC203843L1

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US
Phone: +1-888-267-4436
https://www.origene.com
techsupport@origene.com
EU: info-de@origene.com
CN: techsupport@origene.cn

CALCOCO2 (NM_005831) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: CALCOCO2 (NM 005831) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: CALCOCO2

Synonyms: NDP52

Selection:

Mammalian Cell

None

Vector: pLenti-C-Myc-DDK (PS100064) **E. coli Selection:** Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC203843).

Sequence:

Sgfl-Mlul

Restriction Sites: Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_005831

ORF Size: 1338 bp





CALCOCO2 (NM_005831) Human Tagged Lenti ORF Clone - RC203843L1

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

17q21.32

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 005831.3</u>

 RefSeq Size:
 3682 bp

 RefSeq ORF:
 1341 bp

 Locus ID:
 10241

 UniProt ID:
 Q13137

MW: 52.3 kDa

Cytogenetics:


Gene Summary: This gene encodes a coiled-coil domain-containing protein. The encoded protein functions as

a receptor for ubiquitin-coated bacteria and plays an important role in innate immunity by mediating macroautophagy. Alternatively spliced transcript variants encoding multiple


isoforms have been observed for this gene. [provided by RefSeq, May 2012]



Product images:



Circular map for RC203843L1



Double digestion of RC203843L1 using Sgfl and Mlul