

## Product datasheet for RC203886L1V

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## c20orf72 (MGME1) (NM\_052865) Human Tagged ORF Clone Lentiviral Particle

**Product data:** 

**Product Type:** Lentiviral Particles

**Product Name:** c20orf72 (MGME1) (NM\_052865) Human Tagged ORF Clone Lentiviral Particle

Symbol: c20orf72

Synonyms: bA504H3.4; C20orf72; DDK1; MTDPS11

**Mammalian Cell** 

Selection:

ACCN:

None

**Vector:** pLenti-C-Myc-DDK (PS100064)

Tag: Myc-DDK

ORF Size: 1032 bp

**ORF Nucleotide** 

•

NM 052865

Sequence:

The ORF insert of this clone is exactly the same as(RC203886).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through

naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: <u>NM 052865.2</u>

**RefSeq Size:** 2142 bp

**RefSeq ORF:** 1035 bp

**Locus ID:** 92667

UniProt ID: Q9BQP7

Cytogenetics: 20p11.23

MW: 39.4 kDa







**Gene Summary:** 

The protein encoded by this gene is a nuclear-encoded mitochondrial protein necessary for the maintenance of mitochondrial genome synthesis. The encoded protein is a RecB-type exonuclease and primarily cleaves single-stranded DNA. Defects in this gene have been associated with mitochondrial DNA depletion syndrome-11. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2015]